



Official Website



JUMPER Health

## Shenzhen Jumper Medical Equipment Co., Ltd.

Address: D Building, No. 71, Xintian Road, Fuyong Street, Baoan, Shenzhen, Guangdong, China

Tel: +86 755-26696279 E-mail: info@jumper-medical.com Web: www.jumper-medical.com

## PRODUCT BROCHURE



SINCE  
1998




**JUMPER**

Healthcare within mind

www.jumper-medical.com



# INFRARED THERMOMETER



## ONE KEY MEASUREMENT, ONE SECOND READING.

Infrared Thermometer uses an infrared sensor to detect the infrared radiation emitted by the eardrum and forehead, and converts the value through a microcomputer to obtain the temperature data.

Users can quickly get measurement results after positioning properly the temperature probe in the ear or on the forehead.

In addition to adults' daily temperature measurement, it can also effectively avoid measurement difficulties caused by crying resistance when measuring infants and young children. This kind of thermometer is more suitable for the whole family to use together.

## FEATURES

- Child Mode
- Silent Mode
- Forehead /Ear/Object Mode
- Professional Accuracy
- Hygienic Non-Touch
- 1 Second Reading
- 10s Auto Shut Off
- Smart Fever Alert
- 20 Sets Records



**NEW ARRIVAL**



## JPD-FR302

DISPLAY <b>HD LCD</b>	MEMORY STORAGE <b>20 SETS DATA</b>	TEMPERATURE UNITS <b>°C/°F</b>
AUTOMATIC SHUTDOWN <b>10S ± 1S</b>	BATTERY <b>AAA X 2</b>	ACCURACY (CLINICAL) <b>0.3°C / 0.6°F</b>
MEASURING RANGE		
<b>FOREHEAD MODE</b>	22 - 43 °C / 71.6 - 109.4 °F	
<b>EAR MODE</b>	34 - 43 °C / 93.2 - 109.4 °F	
<b>OBJECT MODE</b>	0 - 100 °C / 32.0 - 212.0 °F	








Forehead/Ear Auto Switch    Magnetic Probe Cap    Child Mode    10s Auto Shut Off    Smart Fever Alert    Silent Mode

**NEW ARRIVAL**



## JPD-FR301

DISPLAY <b>HD LED</b>	MEMORY STORAGE <b>20 SETS DATA</b>	TEMPERATURE UNITS <b>°C/°F</b>
AUTOMATIC SHUTDOWN <b>10S ± 1S</b>	BATTERY <b>AAAX2</b>	ACCURACY (CLINICAL) <b>0.3°C / 0.6°F</b>
MEASURING RANGE		
<b>FOREHEAD MODE</b>	22 - 43 °C / 71.6 - 109.4 °F	
<b>EAR MODE</b>	34 - 43 °C / 93.2 - 109.4 °F	
<b>OBJECT MODE</b>	0 - 100 °C / 32.0 - 212.0 °F	








Forehead/Ear Auto Switch    Magnetic Probe Cap    Child Mode    10s Auto Shut Off    Smart Fever Alert    Silent Mode

**NEW ARRIVAL**



## JPD-FR304

DISPLAY <b>HD LCD</b>	MEMORY STORAGE <b>20 SETS DATA</b>	TEMPERATURE UNITS <b>°C/°F</b>
AUTOMATIC SHUTDOWN <b>10S ± 1S</b>	BATTERY <b>AAAX2</b>	ACCURACY (CLINICAL) <b>0.3°C / 0.6°F</b>
MEASURING RANGE		
<b>EAR MODE</b>	34 - 43 °C / 93.2 - 109.4 °F	








Magnetic Probe Cap    Child Mode    10s Auto Shut Off    Smart Fever Alert    Silent Mode    40 pcs disposable probe cover

**NEW ARRIVAL**



## JPD-FR204

DISPLAY <b>HD LED</b>	MEMORY STORAGE <b>20 SETS DATA</b>	TEMPERATURE UNITS <b>°C/°F</b>
AUTOMATIC SHUTDOWN <b>10S ± 1S</b>	BATTERY <b>AAAX2</b>	ACCURACY (CLINICAL) <b>0.3°C / 0.6°F</b>
MEASURING RANGE		
<b>FOREHEAD MODE</b>	22 - 43 °C / 71.6 - 109.4 °F	
<b>OBJECT MODE</b>	0 - 100 °C / 32.0 - 212.0 °F	







Hygienic Non-Touch    Child Mode    10s Auto Shut Off    Smart Fever Alert    Silent Mode

**NEW ARRIVAL**



## JPD-FR403

DISPLAY <b>HD LCD</b>	MEMORY STORAGE <b>20 SETS DATA</b>	TEMPERATURE UNITS <b>°C/°F</b>
AUTOMATIC SHUTDOWN <b>10S ± 1S</b>	BATTERY <b>AAAX2</b>	ACCURACY (CLINICAL) <b>0.3°C / 0.6°F</b>
MEASURING RANGE		
<b>FOREHEAD MODE</b>	22 - 43 °C / 71.6 - 109.4 °F	
<b>EAR MODE</b>	34 - 43 °C / 93.2 - 109.4 °F	
<b>OBJECT MODE</b>	0 - 100 °C / 32.0 - 212.0 °F	








Forehead/Ear Auto Switch    Magnetic Probe Cap    Child Mode    10s Auto Shut Off    Smart Fever Alert    Silent Mode



## JPD-FR410

DISPLAY <b>HD LCD</b>	MEMORY STORAGE <b>20 SETS DATA</b>	TEMPERATURE UNITS <b>°C/°F</b>
AUTOMATIC SHUTDOWN <b>10S ± 1S</b>	BATTERY <b>AAA X 2</b>	ACCURACY (CLINICAL) <b>0.3°C / 0.6°F</b>
MEASURING RANGE		
<b>FOREHEAD MODE</b>	35 - 42.2 °C / 95 - 108.0 °F	
<b>EAR MODE</b>	35 - 42.0 °C / 95 - 107.6 °F	
<b>OBJECT MODE</b>	0 - 100 °C / 32.0 - 212.0 °F	






Child Mode    10s Auto Shut Off    Smart Fever Alert    Silent Mode



## JPD-FR409-BT

DISPLAY  
**HD LCD**

MEMORY STORAGE  
**20 SETS DATA**

TEMPERATURE UNITS  
**°C/°F**

AUTOMATIC SHUTDOWN  
**10S ± 1S**

BATTERY  
**AAAX2**

ACCURACY (CLINICAL)  
**0.3°C / 0.6°F**

MEASURING RANGE  
**FOREHEAD MODE** 22 - 43 °C / 71.6 - 109.4 °F  
**EAR MODE** 34 - 43 °C / 93.2 - 109.4 °F  
**OBJECT MODE** 0 - 100 °C / 32.0 - 212.0 °F



## JPD-FR300

DISPLAY  
**HD LCD**

MEMORY STORAGE  
**20 SETS DATA**

TEMPERATURE UNITS  
**°C/°F**

AUTOMATIC SHUTDOWN  
**10S ± 1S**

BATTERY  
**AAAX2**

ACCURACY (CLINICAL)  
**0.3°C / 0.6°F**

MEASURING RANGE  
**FOREHEAD MODE** 22 - 43 °C / 71.6 - 109.4 °F  
**EAR MODE** 34 - 43 °C / 93.2 - 109.4 °F  
**OBJECT MODE** 0 - 100 °C / 32.0 - 212.0 °F



## JPD-FR409

DISPLAY  
**HD LCD**

MEMORY STORAGE  
**20 SETS DATA**

TEMPERATURE UNITS  
**°C/°F**

AUTOMATIC SHUTDOWN  
**10S ± 1S**

BATTERY  
**AAAX2**

ACCURACY (CLINICAL)  
**0.3°C / 0.6°F**

MEASURING RANGE  
**FOREHEAD MODE** 35 - 42.2 °C / 95 - 108.0 °F  
**EAR MODE** 35 - 42.0 °C / 95 - 107.6 °F  
**OBJECT MODE** 0 - 100 °C / 32.0 - 212.0 °F



## JPD-FR415

DISPLAY  
**HD LCD**

MEMORY STORAGE  
**20 SETS DATA**

TEMPERATURE UNITS  
**°C/°F**

AUTOMATIC SHUTDOWN  
**10S ± 1S**

BATTERY  
**AAAX2**

ACCURACY (CLINICAL)  
**0.3°C / 0.6°F**

MEASURING RANGE  
**EAR MODE** 34 - 43 °C / 93.2 - 109.4 °F



## JPD-FR400

DISPLAY  
**HD LCD**

MEMORY STORAGE  
**20 SETS DATA**

TEMPERATURE UNITS  
**°C/°F**

AUTOMATIC SHUTDOWN  
**10S ± 1S**

BATTERY  
**AAAX2**

ACCURACY (CLINICAL)  
**0.3°C / 0.6°F**

MEASURING RANGE  
**FOREHEAD MODE** 0 - 100 °C / 32 - 212 °F  
**EAR/OBJECT MODE** 0 - 100 °C / 32 - 212 °F



## JPD-FR401

DISPLAY  
**HD LCD**

MEMORY STORAGE  
**20 SETS DATA**

TEMPERATURE UNITS  
**°C/°F**

AUTOMATIC SHUTDOWN  
**10S ± 1S**

BATTERY  
**AAAX2**

ACCURACY (CLINICAL)  
**0.3°C / 0.6°F**

MEASURING RANGE  
**FOREHEAD MODE** 0 - 100 °C / 32 - 212 °F  
**EAR/OBJECT MODE** 0 - 100 °C / 32 - 212 °F





## JPD-FR100+

DISPLAY  
**HD LCD**

MEMORY STORAGE  
**20 SETS DATA**

TEMPERATURE UNITS  
**°C/°F**

AUTOMATIC SHUTDOWN  
**10S ± 1S**

BATTERY  
**AAAX2**

ACCURACY(CLINICAL)  
**0.3°C / 0.6°F**

MEASURING RANGE  
**FOREHEAD MODE** 32 - 42.2 °C / 89.6 - 107.9 °F  
**EAR/OBJECT MODE** 0.0 - 100 °C / 32.0 - 212.0 °F



## JPD-FR203

DISPLAY  
**HD LCD**

MEMORY STORAGE  
**20 SETS DATA**

TEMPERATURE UNITS  
**°C/°F**

AUTOMATIC SHUTDOWN  
**10S ± 1S**

BATTERY  
**AAAX2**

ACCURACY(CLINICAL)  
**0.3°C / 0.6°F**

MEASURING RANGE  
**FOREHEAD MODE** 35 - 42.2 °C / 95.0 - 108.0 °F  
**OBJECT MODE** 0 - 100.0 °C / 32.0 - 212.0 °F



NEW ARRIVAL



## JPD-FR205

DISPLAY  
**HD LCD**

MEMORY STORAGE  
**20 SETS DATA**


TEMPERATURE UNITS  
**°C/°F**

AUTOMATIC SHUTDOWN  
**10S ± 1S**

BATTERY  
**AAAX2**

ACCURACY(CLINICAL)  
**0.3°C / 0.6°F**

MEASURING RANGE  
**FOREHEAD MODE** 22 - 43 °C / 71.6 - 109.4 °F  
**OBJECT MODE** 0 - 100 °C / 32.0 - 212.0 °F



## JPD-FR202

DISPLAY  
**HD LCD**

MEMORY STORAGE  
**20 SETS DATA**

TEMPERATURE UNITS  
**°C/°F**

AUTOMATIC SHUTDOWN  
**10S ± 1S**

BATTERY  
**AAAX2**

ACCURACY(CLINICAL)  
**0.3°C / 0.6°F**

MEASURING RANGE  
**FOREHEAD MODE** 22 - 43 °C / 71.6 - 109.4 °F  
**OBJECT MODE** 0 - 100 °C / 32.0 - 212.0 °F



NEW ARRIVAL



## FTA10

DISPLAY  
**HD LCD**

MEMORY STORAGE  
**20 SETS DATA**

TEMPERATURE UNITS  
**°C/°F**

AUTOMATIC SHUTDOWN  
**10S ± 1S**

BATTERY  
**AAAX2**

ACCURACY(CLINICAL)  
**0.3°C / 0.6°F**

MEASURING RANGE  
**FOREHEAD MODE** 22 - 43 °C / 71.6 - 109.4 °F  
**OBJECT MODE** 0 - 100 °C / 32.0 - 212.0 °F





Infrared Thermometer															
Parameters															
Product Model	JPD-FR400	JPD-FR409-BT	JPD-FR403	JPD-FR302	JPD-FR301	JPD-FR410		JPD-FR300	JPD-FR401	JPD-FR409	JPD-FR100+	FTA10	JPD-FR205	JPD-FR204	
Display Screen	HD LCD				HD LED	HD LCD		HD LCD					HD LED		
Forehead Mode	0.0-100.0 °C 32.0-212.0 °F	22.0-43.0 °C / 71.6-109.4 °F				35.0-42.0 °C 95.0-108.0 °F		22.0-43.0 °C 71.6-109.4 °F	0.0-100.0 °C 32.0-212.0 °F	35.0-42.0 °C 95.0-108.0 °F	32.0-42.2 °C 89.6-107.9 °F	22.0-43.0 °C / 71.6-109.4 °F			
Ear Mode	0.0-100.0 °C 32.0-212.0 °F	34.0-43.0 °C / 93.2-109.4 °F				35.0-42.0 °C 95.0-107.6 °F		34.0-43.0 °C 93.2-109.4 °F	0.0-100.0 °C 32.0-212.0 °F	35.0-42.0 °C 95.0-107.6 °F	0.0-100.0 °C 32.0-212.0 °F	—			
Object Mode	0.0-100.0 °C / 32.0-212.0 °F							0.0-100.0 °C / 32.0-212.0 °F							
Child Mode	—							—		—					
Silent Mode	—							—		—					
Backlight/Indicator						OPTIONAL: 									—
Clinical Accuracy	± 0.3 °C / ± 0.6 °F							± 0.3 °C / ± 0.6 °F							
Auto Shut Off	10s ± 1s							10s ± 1s							
Features	<ul style="list-style-type: none"> <li>❶ Fever warning</li> <li>❷ °C/°F adjustable</li> <li>❸ Memory storage: 20 sets data</li> <li>❹ Batteries: AAA × 2</li> </ul>		<ul style="list-style-type: none"> <li>❺ Forehead mode: for all ages</li> <li>❻ Ear mode for children: ≥ 3 months old</li> <li>❼ Child mode: ≤ 12 years old</li> </ul>					<ul style="list-style-type: none"> <li>❶ Fever warning</li> <li>❷ °C/°F adjustable</li> <li>❸ Memory storage: 20 sets data</li> <li>❹ Batteries: AAA × 2</li> </ul>		<ul style="list-style-type: none"> <li>❺ Forehead mode: for all ages</li> <li>❻ Ear mode for children: ≥ 3 months old</li> <li>❼ Child mode: ≤ 12 years old</li> </ul>					

## Infrared Thermometer

Parameters



Product Model	JPD-FR202	JPD-FR203	JPD-FR304	JPD-FR415
Display Screen	HD LCD			
Forehead Mode	22.0-43.0 °C 71.6-109.4 °F	35.0-42.2 °C 95.0-108.0 °F	—	
Ear Mode	—	—	34.0-43.0 °C 93.2-109.4 °F	
Object Mode	0.0-100.0 °C 32.0-212.0 °F	—	—	
Child Mode	—	—	—	—
Silent Mode	—	—	—	—
Backlight/Indicator	—	OPTIONAL: 1 ● ● ● ● ● 2 ● ● ● ● ●	—	—
Clinical Accuracy	± 0.3 °C / ± 0.6 °F			
Auto Shut Off	10s ± 1s			
Features	<ul style="list-style-type: none"> <li>1 Fever warning</li> <li>2 °C/°F adjustable</li> <li>3 Memory storage: 20 sets data</li> <li>4 Batteries: AAA × 2</li> <li>5 Forehead mode: for all ages</li> <li>6 Ear mode for children: ≥ 3 months old</li> <li>7 Child mode: ≤ 12 years old</li> </ul>			

