



Official Website



JUMPER Health

## Shenzhen Jumper Medical Equipment Co., Ltd.

**Address:** D Building, No. 71, Xintian Road, Fuyong Street, Baoan, Shenzhen, Guangdong, China

**Tel:** +86 755-26696279 **E-mail:** info@jumper-medical.com **Web:** <https://jumpermed.com>

## PRODUCT BROCHURE

SINCE  
1998









**JUMPER**

Healthcare within mind

<https://jumpermed.com>

# PELVIC FLOOR REHABILITATION DEVICE



-  Bluetooth
-  HD LCD
-  2 separate channels
-  Type-C Charging
-  18 treatment modes
-  Customized mode

## PORTABLE POSTPARTUM REHABILITATION THERAPY DEVICE

BFC10 is a portable Postpartum Rehabilitation Therapy Device, specifically designed for postnatal and menopausal women. It is based on the theory of low frequency electric therapy and pressure biofeedback training, offering a variety of treatments that have been optimized for maximum effectiveness. By stimulate the specific positions of the human body, BFC10 can help with urinary incontinence, repair rectus abdominis muscle, strengthen the vaginal and pelvic floor muscle and provide pelvic floor training to improve the health of women.



NEW  
ARRIVAL



## BFD10

OUTPUT CHANNELS	TREATMENT TIME	POWER SUPPLY
2	5-60 mins	Lithium Battery

NET WEIGHT	DIMENSION
145g	151mm×78mm×32mm

WORKING MODE  
**Continuous operating**

-  HD LCD
-  2 separate channels
-  Type-C Charging
-  18 treatment modes
-  Customized mode

NEW  
ARRIVAL















## BFC10

OUTPUT CHANNELS	TREATMENT TIME	POWER SUPPLY
2	5-60 mins	Lithium Battery

NET WEIGHT	DIMENSION
128g	151mm×82mm×33mm

WORKING MODE  
**Continuous operating**

-  HD LCD
-  2 separate channels
-  Type-C Charging
-  18 treatment modes
-  Customized mode

<div>Pelvic Floor Rehabilitation Device</div> <div>Parameters</div>			
 Product Model	BFD10	BFC10	
 Output Channels	2		
 Treatment Time	5-60 mins		
 Power Supply	Lithium Battery		
 Auto Power Off	✓		
 Low Power Indication	✓		
 Dimension (mm)	151 x 78 x 32	151 x 82 x 33	
 Net Weight (g)	145	128	
 Accessories	Vaginal Electrode Vaginal Pressure Biofeedback Probe Breast Electrode Electrode (100mmx50mm) Electrode (50mmX50mm)	Vaginal Electrode Vaginal Pressure Biofeedback Probe Breast Electrode Electrode (100mmx50mm) Electrode (50mmX50mm)	
 Modes	18		
 Working Mode	Continuous operating		